



S2400 Susie Final Inspection Form
(Ref. "Final Inspection Procedure" PRO-102051-EN-xx)

1. Test Information

1. Simulator Serial No.: _____	5. Type: New ☐ Repair ☐ Upgrade ☐ Show ☐ Loaner ☐
2. Customer: _____	6. Color: Light ☐ Medium ☐ Dark ☐
3. UNI version: _____	7. Power Adapter AC Voltage: _____
4. Options (Special Orders): _____	

2. Cosmetic

Visual Inspection	Pass	Fail	Visual Inspection (continued)	Pass	Fail
1. Skin Condition	☐	☐	8. Interface Panel	☐	☐
2. ECG Inserts	☐	☐	(power button, power connector, RJ45 connector)		
(x4 each of blanks, snaps and patches)			9. Removable Parts:		
3. Defibrillation Inserts	☐	☐	• Ventilation Insert (neck)	☐	☐
(x2 each of blanks and patches)			• Intramuscular (IM) – upper left & right legs	☐	☐
4. Central Line Inserts	☐	☐	• Intramuscular (IM) – left & right buttocks	☐	☐
(x1 each of blank and subclavian central line)			• Left IV forearm	☐	☐
5. Ostomy Inserts – Heal (pink)	☐	☐	• Right IV forearm + Drug Recognition	☐	☐
(x1 each of Ileostomy and Colostomy)			• Glucose Finger (left middle finger)	☐	☐
6. Fill and Drain Ports	☐	☐	• Female Genitalia	☐	☐
(x3 on each side)			• Male Genitalia (Optional)	☐	☐
7. IV Ports	☐	☐	10. Installation of IUD uterus and normal cervix	☐	☐
(x2 on each forearm)			11. Left ulcerated foot	☐	☐

Notes: _____

3. Connection with UNI

Connection Types	Pass	Fail
1. Wired	☐	☐
2. RF Radio	☐	☐
3. Bluetooth	☐	☐

Notes: _____

Use RF Radio connection with UNI for sections 4 – 9.

4. Hardware

Battery	Pass	Fail	Reading	Reservoirs	Pass	Fail
1. Battery status w/o charger	☐	☐	_____	1. Stomach (1.5L)	☐	☐
2. Battery status with charger	☐	☐	_____	2. Gastrostomy (70mL)	☐	☐
				3. Colostomy (1.5L)	☐	☐
				4. Bladder (250mL) + Urinary Catheterization	☐	☐
				5. Central Line (15mL)	☐	☐
				6. CO2 Capability	☐	☐
				7. CO2 Option	☐	☐

Notes: _____



5. Neurological

Seizures	Pass	Fail
1. On	☐	☐
2. Off	☐	☐

Notes: _____

6. Cephalic

Cyanosis	Pass	Fail
1. Intensity 0 – 100%	☐	☐

Notes: _____

7. Ocular

Blink State	Pass	Fail	Pupil Reaction	Pass	Fail
1. Open	☐	☐	1. Left – On/Off	☐	☐
2. Closed	☐	☐	2. Right – On/Off	☐	☐
3. 20 blinks/min	☐	☐			
Pupil Dilation	Pass	Fail			
1. Left (2) / Right (9)	☐	☐			
2. Left (9) / Right (2)	☐	☐			

Notes: _____

8. Airway

Airway	Pass	Fail
1. Tongue Edema – On/Off	☐	☐
2. Laryngospasm – On/Off	☐	☐
3. Pharyngeal Swelling – On/Off	☐	☐

Notes: _____

9. Circulation

Pulses	Pass	Fail	Blood Pressure (using Stethoscope and Blood Pressure cuff)
1. Carotid Right Pulse (right neck)	☐	☐	Systolic (120 +/- 10) Diastolic (80 +/- 10)
2. Carotid Left Pulse (left neck)	☐	☐	Left Arm _____
3. Brachial Right Pulse (inner right bicep)	☐	☐	Right Arm _____
4. Brachial Left Pulse (inner left bicep)	☐	☐	
5. Femoral Right Pulse (right groin)	☐	☐	
6. Femoral Left Pulse (left groin)	☐	☐	
7. Popliteal Right Pulse (back of right knee)	☐	☐	
8. Popliteal Left Pulse (back of left knee)	☐	☐	
9. Radial Pulse (wrist):			
• Left – On/Off	☐	☐	
• Right – On/Off	☐	☐	
10. Pedal Pulse (foot):			
• Left – On/Off	☐	☐	
• Right – On/Off	☐	☐	

Blood Pressure (using Stethoscope and Blood Pressure cuff)
Systolic (120 +/- 10) Diastolic (80 +/- 10)

Left Arm _____
Right Arm _____

Pulse Strength

Pulse Strength Check (Verify ✓) – check either left or right sides

	Carotid	Femoral	Radial	Pedal
BP = 59/40	☐	☐	☐	☐
BP = 69/40	☐	☐	☐	☐
BP = 79/40	☐	☐	☐	☐
BP = 100/40	☐	☐	☐	☐
BP = 200/120	☐	☐	☐	☐
BP = 300/200	☐	☐	☐	☐
BP = 120/80	☐	☐	☐	☐



NIBP	Left Arm	Right Arm	Oxygen Saturation	Pass	Fail
200/120 (+/- 10)	_____	_____	1. O2 Sat. Calibration (left index finger)	☐	☐
160/90 (+/- 10)	_____	_____	2. O2 Sat. Reading (70/80/95 %)	☐	☐
120/80 (+/- 10)	_____	_____			
100/60 (+/- 10)	_____	_____			

Notes: _____

Use Bluetooth connection with UNI for sections 10 – 13.

10. Sounds

Lung Sounds – Anterior (front)	Pass	Fail	Bowel Sounds	Pass	Fail
1. Upper Left Lung	☐	☐	1. Upper Left Bowel	☐	☐
2. Upper Right Lung	☐	☐	2. Upper Right Bowel	☐	☐
3. Lower Left Lung	☐	☐	3. Lower Left Bowel	☐	☐
4. Lower Right Lung	☐	☐	4. Lower Right Bowel	☐	☐
Lung Sounds - Posterior (back)	Pass	Fail	Heart Sounds	Pass	Fail
1. Upper Left Lung	☐	☐	1. Heart Sound (random)	☐	☐
2. Upper Right Lung	☐	☐			
3. Lower Left Lung	☐	☐	Speech	Pass	Fail
4. Lower Right Lung	☐	☐	1. Pre-recorded speech	☐	☐
Throat Sounds	Pass	Fail	2. Streaming Audio with Microphone	☐	☐
1. Throat Sound (random)	☐	☐	3. Streaming Audio with Headphone	☐	☐
			4. Duplex Streaming Audio	☐	☐

Notes: _____

11. Lower Arms / Vascular

Vascular Access	Pass	Fail	Drug Recognition (right forearm)	Pass	Fail
1. Left IV forearm (air test)	☐	☐	1. Identify ID tag	☐	☐
			2. Program drug	☐	☐
			3. Detect 6ml of drug	☐	☐
			4. Delete drug	☐	☐

Notes: _____

12. Breathing

Chest Rise	Pass	Fail	Oral Intubation	Pass	Fail
1. Left – On/Off	☐	☐	1. ET Tube Position 3 (Too Deep) – size #7	☐	☐
2. Right – On/Off	☐	☐	2. Right Mainstem (right lung only)	☐	☐
Respiratory Rate	Pass	Fail	3. ET Tube Position 1 (Normal) – size #7	☐	☐
1. Respiratory Rate (30/20/13 breaths/min)	☐	☐	4. Bilateral chest rise ventilation	☐	☐
			5. ET Tube Position 0 (ET Tube Removed)	☐	☐
			6. Nasal Intubation – size #6	☐	☐
			7. Esophageal Intubation	☐	☐

Notes: _____



13. Cardiac

Cardiac Functions	Pass	Fail	Cardiac Functions (continued)	Pass	Fail
1. Artificial Ventilation with BVM mask	☐	☐	6. Rhythms and Heart Rates OK	☐	☐
2. Auto Lung Sound	☐	☐	7. ECG Rhythm at Shock Sites	☐	☐
3. Chest Compression	☐	☐	8. Sternum/Apex Pacing Detection	☐	☐
4. Auto Carotid Pulse	☐	☐	9. Sternum/Apex Pacing Capture	☐	☐
5. ECG Rhythm at ECG Sites (4-lead)			10. Defibrillation with snap cable (150J)	☐	☐
• With snaps	☐	☐	11. Install blank ECG and Defibrillation inserts	☐	☐
• With patches	☐	☐			

Notes: _____

Disconnect from UNI for section 14.

14. External Features

Articulation (manual movement)	Pass	Fail	Waist Rod Knobs	Pass	Fail
1. Left shoulder	☐	☐	1. Check left and right knobs are fastened onto the waist rod	☐	☐
2. Right shoulder	☐	☐			
3. Left hip	☐	☐	Labels	Pass	Fail
4. Right hip	☐	☐	1. RF Module (USB cable)	☐	☐
5. Neck	☐	☐	2. Tablet/Desktop (lower left corner)	☐	☐
Remove and Attach Lower Legs	Pass	Fail	3. Charger (AC power adapter)	☐	☐
1. Right leg	☐	☐	4. S2400 serial number (posterior insert)	☐	☐
2. Left leg	☐	☐			

Notes: _____

15. Observations

16. Inspection Information

Inspected by: _____ Date: _____

----- Document Information -----

Document No.: FORM-102239-EN-01	Edited by: Perry Tai	Date: 2/22/24	Change(s): Initial released version.
------------------------------------	-------------------------	------------------	---